






MEDICAL SONITA MEDIC SURGICAL SUPPLIES

OT Consumables Division



Gynaecology Solutions



| www.sonitamedic.com.au |

■ EMAIL: customers@sonitamedic.com.au

■ PHONE: 1800 SONITA (766 482)

■ MOBILE: 0423 178 332

SYDNEY :

Suite 242 / 60
Walker St, North Sydney,
NSW 2060

MELBOURNE :

Suite 121 / 358 Clarendon St,
South Melbourne,
VIC 3205

BRISBANE :

524 Old Bay Road
Burrpengary East,
QLD 4505

PERTH :

Suite 76 Ground Floor
50, St Georges Terrace,
Perth WA 6000

STERILE SINGLE USE VAGINAL SPECULUM WITH SMOKE EXTRACTOR




- Sterile Single Use Vaginal Speculum with smoke evacuation tube.
- Designed for use with Leetz loop diathermy.
- Smooth round edges increase patient comfort.
- Colour efficiently reflects light, ensuring maximum visibility for the procedure.
- Improved configuration of holes which removes smoke effectively.
- Proven strength, with Failsafe Points uniquely built into each vaginal speculum.
- Supplied sterile in box of 25.



◆ broad size

◆ medium long size

◆ medium size



Code	Description	Size
• PEL-400111	Vaginal Speculum with Smoke Extractor	Medium
• PEL-400112	Vaginal Speculum with Smoke Extractor	Medium/Long
• PEL-400113	Vaginal Speculum with Smoke Extractor	Broad

STERILE DISPOSABLE SPACKMAN UTERINE CANNULA



- Provides a cost-effective way of guaranteeing sterility and preventing any chance of cross-contamination.
- Manufactured from Stainless Steel and Lexan Polymer.
- Good seal on cervix.
- Simple locking mechanism holds fast and it securely holds the vellsullum forceps in place.
- Luer connection on proximal end for fitting to syringe
- Latex Free
- Single Use
- Supplied Sterile in box of 25



The unique locking device holds the vellsullum forceps in place.

Acorn to achieve an effective seal.



Code

Description


• **SHEF-UC001** Sterile Single Use Uterine Cannula

LAPAROSCOPIC CONSUMABLES


STERILE SINGLE-USE PTFE COATED NON-STICK




- Minimize adhesion & prevent bleeders from resuming.
- Safer for the patient & reduces the risk of damaging surrounding tissue.
- Electrode does not need to be removed and cleaned.
- Latex free & supplied in box of 10.



● Spatula Tip




● Needle Tip




● Hook Tip

Code	Description	Size
● BM-LES184	Spatula Electrode, with Spatula Diathermy	5mm/330mm long
● BM-LES224	Needle Electrode	5mm/330mm long
● BM-LES104	Hook Electrode, with L-Wire Diathermy	5mm/330mm long
● BM-LEL104	Hook Electrode, with L-Wire Diathermy	5mm/450mm long

STERILE SINGLE USE L HOOK WITH CABLE



● Spatula Tip



● Needle Tip



● Hook Tip

Code	Description	Size
● BM-HE002	L Hook Electrode	5mm/330mm long
● BM-HE002-C	L Hook Electrode, with 3m cable	5mm/330mm long
● BM-SE002	Spatula Electrode	5mm/330mm long
● BM-SE002-C	Spatula Electrode, with 3m cable	5mm/330mm long
● BM-HE007	L Hook Electrode, monopolar diathermy connection	5mm/330mm long
● BM-HE007-C	L Hook Electrode, monopolar diathermy connection, with 3m cable	5mm/330mm long
● BM-SE007-C	Spatula Electrode, monopolar diathermy connection, with 3m cable	5mm/330mm long

MEDICAL

■ PHONE:

1800 SONITA (766 482)

■ MOBILE:

0423 178 332

SONITA MEDIC SURGICAL

■ EMAIL: customers@sonitamedic.com.au

■ WEB: www.sonitamedic.com.au

SUPPLIES

■ OFFICE's in:

• SYDNEY


• MELBOURNE

• BRISBANE

• PERTH

STERILE SINGLE USE SMOKE EVACUATION PENCILS

STERILE SINGLE-USE SMOKE EVACUATION PENCILS



- Smoke Pencils Extension Suction Tip
70/125/165/200mm long

- Offers an all-in-one design to eliminate the need for multiple smoke evacuation products that need to be assembled in the sterile field.
- Light flexible tubing allows axial rotation.
- Designed for maximum visibility in the surgical site and comfort in the surgeon's hand.
- Compatible with a wide range of electrodes.
- Greater capture of harmful pathogens contained in surgical smoke.

Code	Description
● MED-04-301	Smoke Evacuation Pencil, Push Button, Handcontrol, 3 Prong Switch, Removable PTFE non-stick coated 70mm Blade Electrode, 3m Cable, Sterile Single Use
● MED-04-302	Smoke Evacuation Pencil, Push Button, Handcontrol, 3 Prong Switch, Removable Stainless Steel 70mm Blade Electrode, 3m Cable, Sterile Single Use
● MED-04-303	Smoke Evacuation Pencil, Push Button, Handcontrol, Cradle and Tubing, 3 Prong Switch, Removable PTFE non-stick coated 70mm Blade Electrode, 3m Cable, Sterile Single Use
● MED-04-304	Smoke Evacuation Pencil, Push Button, Handcontrol, 3 Prong Switch, Cradle and Tubing, Removable Stainless Steel 70mm Blade Electrode, 3m Cable, Sterile Single Use

** Add **-H** to code if you require Smoke Evacuation Pencil with Holster
Add **-R** to code if you require Smoke Evacuation Pencil with Rocker Switch
Add **-RH** to code if you require Smoke Evacuation Pencil with both Holster and Rocker Switch

Code	Description
● MED-04-305	Smoke Pencils Extension Suction Tip, 70mm long, Sterile Single Use
● MED-04-306	Smoke Pencils Extension Suction Tip, 125mm long, Sterile Single Use
● MED-04-307	Smoke Pencils Extension Suction Tip, 165mm long, Sterile Single Use
● MED-04-308	Smoke Pencils Extension Suction Tip, 200mm long, Sterile Single Use
● MED-04-309	Smoke Evacuation Tubing, 3/8 inch x 10 Feet
● MED-04-309A	Smoke Evacuation Tubing, 7/8 inch x 10 Feet

COTSWOLD MICRO DISSECTION NEEDLES


Designed for use with virtually all electrosurgery generators and the majority of electrosurgery handles / finger switches, the Cotswold range of Monopolar micro-dissection needles are exceptionally fine needles for precise Monopolar dissection and offers the surgeon the ability to direct the power from the electrosurgery generator to an area (typically 3 micron) at the needle tip.

Traditionally produced "sharp" electrosurgery needles become blunt very quickly during a procedure due to the stainless steel "melting", particularly with such a fine point. The Cotswold needles are produced in rigid, high temperature pure tungsten which has a melt temperature of 3400°C, significantly higher than those temperatures used in electrosurgery.

BM-1650




BM-1651



BM-1654



BM-1692



BM-1693



BM-1694



Code	Description
• BM-1650	COTSWOLD Micro Dissection Needle, Single Use, Sharp Tungsten Tip, 30mm Shaft, 10mm Tip
• BM-1651	COTSWOLD Micro Dissection Needle, Single Use, Sharp Tungsten Tip, 30mm Shaft, 20mm Tip
• BM-1654	COTSWOLD Micro Dissection Needle, Single Use, Sharp Tungsten Tip, 55mm Shaft, 10mm Tip
• BM-1692	COTSWOLD Micro Dissection Needle, Single Use, Sharp Tungsten Tip, 30mm Shaft, 5mm Angled Tip
• BM-1693	COTSWOLD Micro Dissection Needle, Single Use, Sharp Tungsten Tip, 40mm Shaft, 5mm Angled Tip
• BM-1694	COTSWOLD Micro Dissection Needle, Single Use, Sharp Tungsten Tip, 55mm Shaft, 5mm Angled Tip

BLADE ELECTRODES

- Non-Stick PTFE coating minimizes adhesion and prevents bleeders from resuming when concluding the procedure.
- This allows for fulguration to be achieved with greater success rates than with standard stainless steel ball electrodes.
- Single-use and supplied sterile in box of 10.

BM-2950



BM-2920



BM-2951



BM-2921



BM-2926



BM-2925



BM-2928



BM-2927



BM-2933



BM-2929



Code	Description	Size
• BM-2950	Blade Electrode Shielded Tip	64mm
• BM-2920	Blade Electrode	64mm
• BM-2951	Blade Electrode Shielded Tip	70mm
• BM-2921	Blade Electrode	70mm
• BM-2926	Blade Electrode Shielded Tip	100mm

Code	Description	Size
• BM-2925	Blade Electrode	100mm
• BM-2928	Blade Electrode Shielded Tip	130mm
• BM-2927	Blade Electrode	130mm
• BM-2933	Blade Electrode Shielded Tip	165mm
• BM-2929	Blade Electrode	165mm

STERILE SINGLE USE BALL/NEEDLE ELECTRODES

SINGLE-USE 2.4MM SHAFT



■ **BM-16-010-07**
Blade Electrodes 70mm



■ **BM-16-012-07**
Blade Electrodes (cur) 70mm



■ **BM-16-020-07**
Needle Electrodes 70mm



■ **BM-16-022-07**
Needle Electrodes (cur) 70mm



Code	Description	Length
STRAIGHT		
● BM-16-030-02	2mm Ball Electrode	70mm
● BM-16-030-03	3mm Ball Electrode	70mm
● BM-16-030-04	4mm Ball Electrode	70mm
● BM-16-030-05	5mm Ball Electrode	70mm
● BM-16-030-06	6mm Ball Electrode	70mm
● BM-16-030-07	7mm Ball Electrode	70mm
● BM-16-030-08	8mm Ball Electrode	70mm
CURVED		
● BM-16-032-02	2mm Ball Electrode	70mm
● BM-16-032-03	3mm Ball Electrode	70mm
● BM-16-032-04	4mm Ball Electrode	70mm
● BM-16-032-05	5mm Ball Electrode	70mm
● BM-16-032-06	6mm Ball Electrode	70mm
● BM-16-032-07	7mm Ball Electrode	70mm
● BM-16-032-08	8mm Ball Electrode	70mm

STERILE SINGLE USE BALL/NEEDLE ELECTRODES

SINGLE-USE 2.4MM SHAFT



■ **BM-16-040-13**

Blade Electrodes 130mm

■ **BM-16-042-13**

Blade Electrodes (Cur) 130mm



■ **BM-16-050-13**

Needle Electrodes 130mm

■ **BM-16-052-13**

Needle Electrodes (Cur) 130mm



Code	Description	Length
STRAIGHT		
● BM-16-060-02	2mm Ball Electrode	130mm
● BM-16-060-03	3mm Ball Electrode	130mm
● BM-16-060-04	4mm Ball Electrode	130mm
● BM-16-060-05	5mm Ball Electrode	130mm
● BM-16-060-06	6mm Ball Electrode	130mm
● BM-16-060-07	7mm Ball Electrode	130mm
● BM-16-060-08	8mm Ball Electrode	130mm

SINGLE-USE 4.0MM SHAFT




■ **BM-17-110-07**

Blade Electrode 4.0mm ø




■ **BM-17-111-07**

Blade Electrodes (Cur) 4.0mm ø



■ **BM-17-120-07**

Needle Electrodes 4.0mm ø



■ **BM-17-121-07**

Needle Electrodes (Cur) 4.0mm ø

STERILE SINGLE USE BALL/NEEDLE ELECTRODES

SINGLE-USE 4.0MM SHAFT



Code	Description	Length
STRAIGHT		
• BM-17-130-02	2mm Ball Electrode	70mm
• BM-17-130-03	3mm Ball Electrode	70mm
• BM-17-130-04	4mm Ball Electrode	70mm
• BM-17-130-05	5mm Ball Electrode	70mm
• BM-17-130-06	6mm Ball Electrode	70mm
• BM-17-130-07	7mm Ball Electrode	70mm
• BM-17-130-08	8mm Ball Electrode	70mm
CURVED		
• BM-17-131-02	2mm Ball Electrode	70mm
• BM-17-131-03	3mm Ball Electrode	70mm
• BM-17-131-04	4mm Ball Electrode	70mm
• BM-17-131-05	5mm Ball Electrode	70mm
• BM-17-131-06	6mm Ball Electrode	70mm
• BM-17-131-07	7mm Ball Electrode	70mm
• BM-17-131-08	8mm Ball Electrode	70mm

SINGLE-USE 4.0MM SHAFT



- **BM-17-140-13**
Blade Electrodes 4.0mm ø
130mm
- **BM-17-142-13**
Blade Electrodes (Cur) 4.0mm ø
130mm



- **BM-17-150-13**
Needle Electrodes 4.0mm ø
130mm
- **BM-17-152-13**
Needle Electrodes (Cur) 4.0mm ø
130mm

Code	Description	Length
STRAIGHT		
• BM-17-160-02	2mm Ball Electrode 4.0mm ø	130mm
• BM-17-160-03	3mm Ball Electrode 4.0mm ø	130mm
• BM-17-160-04	4mm Ball Electrode 4.0mm ø	130mm
• BM-17-160-05	5mm Ball Electrode 4.0mm ø	130mm
• BM-17-160-06	6mm Ball Electrode 4.0mm ø	130mm
• BM-17-160-07	7mm Ball Electrode 4.0mm ø	130mm
• BM-17-160-08	8mm Ball Electrode 4.0mm ø	130mm


STERILE SINGLE USE BALL/NEEDLE ELECTRODES

SINGLE-USE 1.6MM SHAFT



■ BM-18-210-01

Blade Insulated Electrodes 70mm long




■ BM-18-210-02

Needle Electrodes (str) 70mm long




■ BM-18-210-03

Needle Electrodes (cur) 70mm long




■ BM-18-210-04

Ball Electrodes 70mm long



■ BM-18-210-05

Diamond Loop Electrodes 70mm long




■ BM-18-210-06

Circle Loop Electrodes 70mm long


STERILE SINGLE USE LOOP ELECTRODES

SINGLE-USE 2.4MM SHAFT




T-BAR LOOP ELECTRODE

Code	Description	Length
• BM-19-310-01T	T Bar Loop Electrode 20 x 8mm shaft	130mm
• BM-19-310-02T	T Bar Loop Electrode 20 x 10mm shaft	130mm
• BM-19-310-03T	T Bar Loop Electrode 20 x 15mm shaft	130mm
• BM-19-310-04T	T Bar Loop Electrode 20 x 18mm shaft	130mm
• BM-19-310-05T	T Bar Loop Electrode 20 x 20mm shaft	130mm
• BM-19-310-06T	T Bar Loop Electrode 25 x 10mm shaft	130mm
• BM-19-310-07T	T Bar Loop Electrode 25 x 15mm shaft	130mm
• BM-19-310-08T	T Bar Loop Electrode 25 x 20mm shaft	130mm



SAIL FRAME LOOP ELECTRODE

Code	Description	Length
• BM-19-320-01SA	Sail Frame Loop Electrode 2.4 x 10mm shaft	130mm
• BM-19-320-02SA	Sail Frame Loop Electrode 2.4 x 15mm shaft	130mm
• BM-19-320-03SA	Sail Frame Loop Electrode 2.4 x 18mm shaft	130mm
• BM-19-320-04SA	Sail Frame Loop Electrode 2.4 x 20mm shaft	130mm




SQUARE LOOP ELECTRODE

Code	Description	Length
• BM-19-350-01SQ	Square Loop Electrode 10mm shaft	130mm
• BM-19-350-02SQ	Square Loop Electrode 15mm shaft	130mm
• BM-19-350-03SQ	Square Loop Electrode 20mm shaft	130mm
• BM-19-350-04SQ	Square Loop Electrode 25mm shaft	130mm


STERILE SINGLE USE LOOP ELECTRODES

SINGLE-USE 2.4MM SHAFT




STANDARD LOOP ELECTRODE

Code	Description	Length
• BM-19-330-01L	Standard Loop Electrode 8mm shaft	130mm
• BM-19-330-02L	Standard Loop Electrode 10mm shaft	130mm
• BM-19-330-03L	Standard Loop Electrode 12mm shaft	130mm
• BM-19-330-04L	Standard Loop Electrode 15mm shaft	130mm
• BM-19-330-05L	Standard Loop Electrode 18mm shaft	130mm
• BM-19-330-06L	Standard Loop Electrode 20mm shaft	130mm
• BM-19-330-07L	Standard Loop Electrode 22mm shaft	130mm
• BM-19-330-08L	Standard Loop Electrode 25mm shaft	130mm



LOW-PROFILE LOOP ELECTRODE

Code	Description	Length
• BM-19-340-01LP	Low-profile Loop Electrode 8mm shaft	130mm
• BM-19-340-02LP	Low-profile Loop Electrode 10mm shaft	130mm
• BM-19-340-03LP	Low-profile Loop Electrode 12mm shaft	130mm
• BM-19-340-04LP	Low-profile Loop Electrode 15mm shaft	130mm
• BM-19-340-05LP	Low-profile Loop Electrode 18mm shaft	130mm
• BM-19-340-06LP	Low-profile Loop Electrode 20mm shaft	130mm
• BM-19-340-07LP	Low-profile Loop Electrode 22mm shaft	130mm
• BM-19-340-08LP	Low-profile Loop Electrode 25mm shaft	130mm



TRAFFORD LOOP ELECTRODE

Code	Description	Length
• BM-19-360-01TR	Trafford Loop Electrode 10mm shaft	130mm
• BM-19-360-02TR	Trafford Loop Electrode 15mm shaft	130mm
• BM-19-360-03TR	Trafford Loop Electrode 20mm shaft	130mm
• BM-19-360-04TR	Trafford Loop Electrode 25mm shaft	130mm

LLETZ LOOPS ELECTRODES

- Ultra-fine rigid tungsten wire for strength and minimum eschar build-up.
- Protective blister packaging.
- Supplied sterile in box of 5.
- Non-Slip Shaft grip for ease of removal.



BM-1740



BM-1750



BM-1742



BM-1746



BM-1743



BM-1751

Code	Description	Size
● BM-1740	Loop Electrode Pink	10mm x 5mm Loop
● BM-1750	Loop Electrode Aqua	10mm x 8mm Loop
● BM-1742	Loop Electrode Green	10mm x 10mm Loop
● BM-1746	Loop Electrode Orange	15mm x 8mm Loop
● BM-1743	Loop Electrode Blue	15mm x 15mm Loop
● BM-1751	Loop Electrode Grey	20mm x 12mm Loop

LLETZ LOOPS ELECTRODES


- Ultra-fine rigid tungsten wire for strength and minimum eschar build-up.
- Protective blister packaging.
- Supplied sterile in box of 5.
- Non-Slip Shaft grip for ease of removal.



BM-1752



BM-1745



BM-1744



BM-1735



BM-1731



BM-1732

Code	Description	Size
● BM-1752	Loop Electrode Yellow	20mm x 13mm Loop
● BM-1745	Loop Electrode Red	20mm x 15mm Loop
● BM-1744	Loop Electrode Beige	20mm x 20mm Loop
● BM-1735	Square Loop Electrode Light Green	5mm x 5mm Loop
● BM-1731	Square Loop Electrode Green	10mm x 4mm Loop
● BM-1732	Square Loop Electrode Light Brown	10mm x 10mm Loop

MEDICAL SONITA MEDIC SURGICAL SUPPLIES

| www.sonitamedic.com.au |

■ EMAIL: customers@sonitamedic.com.au

■ PHONE: 1800 SONITA (766 482)

■ MOBILE: 0423 178 332

SYDNEY :

Suite 242 / 60
Walker St, North Sydney,
NSW 2060

MELBOURNE :

Suite 121 / 358 Clarendon St,
South Melbourne,
VIC 3205

BRISBANE :

524 Old Bay Road
Burrpengary East,
QLD 4505

PERTH

Suite 76 Ground Floor
50, St Georges Terrace,
Perth WA 6000