SONITA MEDIC SURGICAL

SUPPLIES

STERILE SINGLE USE SPACKMAN UTERINE CANNULA





Designed to replace reusable 'Spackman' style cannula, comes with a Luer fitting, can be used to introduce dye and other fluids to the uterus, and is often needed in the treatment of tubal occlusion. Simple yet effective.

Disposable Spackman Uterine Cannula

- Incorporates a unique locking device for ease of use.
- Allows the forcep to be positioned on the cannula without using the finger rings.
- Materials: Stainless Steel and Lexan Polymer
- Sterilisation: ETO
- Comformity: CEO120
- Latex: None
- Individual packed in double peel pouch
- Box of 25

MEDICAL

■ PHONE:

1800 SONITA (766 482)

MOBILE:

0423 178 332

SONITA MEDIC SURGICAL

■ EMAIL: customers@sonitamedic.com.au

■ WEB: www.sonitamedic.com.au

SUPPLIES

OFFICE's in:

. SYDNEY

. MELBOURNE

. BRISBANE

. PERTH