


SILICONE MATS

BECAUSE DELICATE ITEMS NEED PROTECTION

MULTIPAK®

FINGER MATS: The Riley Medical Finger Mats provide a convenient method of positioning and protecting surgical instruments during storage, transportation, and sterilization. Four finger configurations allow utilization of the mats with a variety of instruments.



Underside of small Diamond Mat

FEATURES

- Silicone Rubber
- 4 Styles; Multiple Sizes
- Fine Pin
- Unique Diamond Grid Pattern (2 size grids) Large & Small
- Small Diamond Grid features raised underside
- Coarse Pin
- Large Sheet

BENEFITS

- Withstands repeated autoclaving
- Protects all types of instruments
- Closely spaced flexible pins adapt to microsurgical instruments
- Superior steam penetration & fast drying
- Makes cleaning & drying a snap
- Medium spaced semirigid pins are suitable for less delicate instruments
- Large sheet allows the mat to be cut to fit existing tray or baskets

WA Distributor:



**RILEY
MEDICAL**

*Manufacturing a full line
of engineered trays
that protect and organize
your instruments*

FLASHPAK®
*Mobile Sterilization
Container System*

MICROPAK®
Microsurgical Instrument Trays


MULTIPAK®
Instrument Organizing Trays

SCOPEPAK®
Endoscope Trays

TRANSPAK®
*Instrument Soaking
& Transport Trays*

MULTIPAK®

SILICONE MATS



choose a style ➤

For MultiPak tray:	Size inches (mm)
5000 series	8.9 x 19.7" (226 x 500)
5000 series w/ wire rack	8.9 x 18" (226 x 460)
5010 series	8.4 x 19.2" (213 x 487)
5012 series	3.7 x 19.2" (94 x 487)
5022 series	9.1 x 8.4" (231 x 213)
5023 series	8.4 x 12.5" (213 x 317)
5006 series	5.4 x 12.5" (137 x 317)

Fine Pin	Sm. Diamond	Lg. Diamond	Coarse Pin
Cat.#	Cat.#	Cat.#	Cat.#
5000-04	5000-04SD	5000-04D	30P*
5110-04	5110-04SD	5110-04D	30P*
5010-04	5010-04SD	5010-04D	30P*
5012-04	5012-04SD	5012-04D	30P*
5022-04	5022-04SD	5022-04D	30P*
5023-04	5023-04SD	5023-04D	30P*
5006-04	NA	5006-04D	30P*

*Some trimming of this mat may be required, measures 10 x 18" (254 x 460mm), also available in a fine pin and unperforated style(s).

To order mats for Riley trays not listed, trays using wire instrument racks or specialty trays please consult your Riley Distributor for details.

WA Distributor: